

M2000 SERIES PLATINUM SILICONES



Custom made using only premium quality platinum-cure silicone materials, M2000 earmolds are excellent choice for most difficult cases.

Recommended Applications: BTE Hearing Aids / Power BTE Hearing Aids / RIC/RITE
Pediatrics / Adults / Senior

- Comfort
- Hypoallergenic
- Excellent resilience
- Minimum discoloration
- High longevity
- High flexibility
- Perfect acoustic seal
- Thousands of color combinations

Material	M35 (formerly known as M2000)	M25 NEW	M45 NEW
Hardness	Shore A 37	Shore A 23	Shore A 45
Recommended age group	Pediatrics / Adults	Pediatrics / Adults	Senior
Hypoallergenic	YES	YES	YES
Default color	Colorless Clear	Colorless Clear	Colorless Clear



For further information, please call us at **1-800-523-7672** or send us an email at **info@microsonic-inc.com** or visit online at **www.microsonic-inc.com**
 ©2017 Microsonic Inc. All Rights Reserved. Woman owned and operated. ISO 9001:2008 Certified

