

M2000 Series

Platinum Silicones



Exceptional comfort for all day use.

Custom made using only premium quality platinum-cure silicone materials, M2000 earmolds are an excellent choice for even the most difficult cases.

- Comfortable
- Excellent resilience
- Hypoallergenic
- Minimum discolorization
- Product longevity
- Highly flexible
- Perfect acoustic seal
- Thousands of color combinations



www.microsonic-inc.com
 Microsonic nc,
 2960 Duss Avenue,
 Ambridge, PA 15003
 P: 855-840-6231
 F: 724-266-6309
 info@microsonic-inc.com



	M25	M35 (formerly M2000)	M45
Hardness	Shore A 23	Shore A 37	Shore A 45
Recommended age group	Pediatrics / Adults	Pediatrics / Adults	Senior
Hypoallergenic	YES	YES	YES
Default color	Clear	Clear	Clear

For more information, please call us at 1-800-523-7672
 Email: info@microsonic-inc.com or Visit: microsonic-inc.com
 Reserved. Woman owned and operated. ISO 9001 Certified
 ©2018 Microsonic Inc. All Rights Reserved.