

# M2000 Series

## Platinum Silicones



### Exceptional comfort for all day use.

Custom made using only premium quality platinum-cure silicone materials, M2000 earmolds are an excellent choice for even the most difficult cases.

- Comfortable
- Excellent resilience
- Hypoallergenic /Antibacterial
- Minimum discolorization
- Product longevity
- Highly flexible
- Perfect acoustic seal
- Thousands of color combinations



www.microsonic-inc.com  
 Microsonic inc.  
 2960 Duss Avenue,  
 Ambridge, PA 15003  
 P: 855-840-6231  
 F: 724-266-6309  
 info@microsonic-inc.com



	M25	M35 (formerly M2000)	SlikFit®	M45
Hardness	Shore A 23	Shore A 37	Shore A 37	Shore A 45
Recommended age group	Pediatrics / Adults	Pediatrics / Adults	All Ages	Senior
Hypoallergenic	YES	YES	YES	YES
Antibacterial	YES	YES	YES	YES
Default color	Clear	Clear	Translucent	Clear
Surface Finish	Glossy	Glossy	Glossy	Silky / Matte

For more information, please call us at 1-800-523-7672  
 Email: [info@microsonic-inc.com](mailto:info@microsonic-inc.com) or Visit: [microsonic-inc.com](http://microsonic-inc.com)  
 Reserved. Woman owned and operated. ISO 9001 Certified 2015  
 ©2018 Microsonic Inc. All Rights Reserved.