


## Hearing Protection PEDIATRIC

### SoftGuard™ Super Soft Silicone



Mean Attenuation Values for SoftGuard™ Custom Hearing Protection Earmolds



#### Exceptional comfort for all day use.

- NRR 26, Mean 40 dB, Attenuation up to 48 dB
- Ideal for active children / Maintains acoustic seal
- Permits hearing conversation
- Super soft for maximum comfort / Hardness Shore A 30
- Hypoallergenic & Antibacterial
- Thousands of color combinations to choose from
- 3D scanned & archived: Hassle-free re-orders/remakes
- Available Styles: Shell, Canal, and Canal-Shell (1/2 Shell)

www.microsonic-inc.com  
Microsonic nc.  
2960 Duss Avenue,  
Ambridge, PA 15003  
P: 855-840-6231  
F: 724-266-6309  
info@microsonic-inc.com

